

VEQ 360°



Una risorsa indispensabile per confermare e migliorare la qualità di tutto il tuo laboratorio.

Catalogo 2026



INSTAND

Liquido Cefalo Rachidiano 2026

Cerebrospinal Fluid Analysis - CSF Ferritin

Uhr / Zerr

467

Samples : 1 CSF sample of 1 ml
Sample properties: CSF
Shipping Information: Shipping at ambient temperature
Certificate Validity: 24 months

Ferritin - CSF

Neurochemical diagnostics of hemorrhages

May (3)

Registration deadline:	27.02.26
Sample shipment:	28.04.26
Closing Date:	15.05.26

Cerebrospinal Fluid Analysis - Dementia Markers

Uhr / Zerr

466

Samples : 1 sample of 1 ml Sample properties: CSF Shipping Information: Shipping at ambient temperature Certificate Validity: 12 months		
Amyloid β Peptid 40 (Aβ40) Amyloid β Peptid 42 (Aβ42) Amyloid β Peptid Ratio (Aβ 42/40) PhosphoTau Protein (pTau) Total Tau Protein (Tau) Neurochemical dementia diagnostics: Support for differential diagnosis of dementia by interpretation of the combined analysis of different marker proteins in CSF.		
	May (3)	Oct. (6)
Registration deadline:	27.02.26	07.08.26
Sample shipment:	28.04.26	06.10.26
Closing Date:	15.05.26	23.10.26

Cerebrospinal Fluid Analysis - Free light chains Kappa

Uhr / Zerr

468

Samples : 2 sets = CSF + serum of 1 ml Sample properties: CSF + serum Shipping Information: Shipping at ambient temperature Certificate Validity: 24 months	
kappa-free light chains - CSF kappa-free light chains - serum Neurochemical diagnostics of inflammatory CNS processes	
	May (3)
Registration deadline:	27.02.26
Sample shipment:	28.04.26
Closing Date:	15.05.26

Cerebrospinal Fluid Analysis - Lactate/ Glucose

Uhr / Zimmermann

464

Samples : 2 samples of 0.5 ml Sample properties: CSF + Serum (lyophilized) Shipping Information: Shipping at ambient temperature Certificate Validity: 12 months / 6 months				
Glucose - CSF Glucose quotient (CSF/Serum) Glucose quotient - diagnostic comment Interpretation of the report Lactate Lactate - diagnostic comment Calculation of the CSF/serum ratio of glucose. Interpretation with respect to the reference range.				
	Jan. (1)	May (3)	Aug. (5)	Oct. (6)
Registration deadline:	14.11.25	27.02.26	12.06.26	07.08.26
Sample shipment:	13.01.26	28.04.26	11.08.26	06.10.26
Closing Date:	30.01.26	15.05.26	28.08.26	23.10.26

participation according to Rili-BÄK:
quarterly

Samples : 4 samples of 1ml
 Sample properties: CSF + Serum
 Shipping Information: Shipping at ambient temperature
 Certificate Validity: 12 months

Antibodies against Measles viruses, Rubella viruses, Varicella zoster viruses, Herpes simplex viruses (1+2).
 Analysis of virus-specific antibodies in CSF and serum for detection of specific antibody response in central nervous system (MRZ reaction).

Detection of titers is not suitable. Evaluation as antibody index and interpretation of combined data.

	May (3)	Oct. (6)
Registration deadline:	27.02.26	07.08.26
Sample shipment:	28.04.26	06.10.26
Closing Date:	15.05.26	23.10.26

Samples : 2 samples of 0.5ml
 Sample properties: CSF + Serum
 Shipping Information: Shipping at ambient temperature
 Certificate Validity: 12 months

Borrelia IgG Antibody Index (AI) calculation Interpretation of the report

Detection of titers is not suitable. Calculation of antibody index (Certificate) and interpretation (Certificate of Participation).

	May (3)	Oct. (6)
Registration deadline:	27.02.26	07.08.26
Sample shipment:	28.04.26	06.10.26
Closing Date:	15.05.26	23.10.26

Samples : 2 samples of 150 µl
 Sample properties: CSF + Serum
 Shipping Information: Shipping at ambient temperature
 Certificate Validity: 12 months

participation according to Rili-BÄK:
 bi-annual

Interpretation of the report Oligoclonal IgG

The acquisition of samples is demanding, as it is not possible to use pooled CSF, i.e. the samples are from individual patients (catheter, drainage, etc.).

	May (3)	Oct. (6)
Registration deadline:	27.02.26	07.08.26
Sample shipment:	28.04.26	06.10.26
Closing Date:	15.05.26	23.10.26

Samples : 4 samples of 1ml
 Sample properties: CSF + Serum
 Shipping Information: Shipping at ambient temperature
 Certificate Validity: 6 months

**participation according to Rili-BÄK:
 quarterly**

(Evaluation of CSF/serum quotients in diagrams and age-related interpretation)

Albumin
 IgG - CSF

Albumin - CSF
 IgM

IgA IgA - CSF
 IgM - CSF

IgG
 Protein - CSF

Sample volumes in the CSF proficiency tests (No. 460-466) are limited as it is difficult to acquire large volumes of CSF.

	Jan. (1)	May (3)	Aug. (5)	Oct. (6)
Registration deadline:	14.11.25	27.02.26	12.06.26	07.08.26
Sample shipment:	13.01.26	28.04.26	11.08.26	06.10.26
Closing Date:	30.01.26	15.05.26	28.08.26	23.10.26